

MacOS compatibility chart DinoXcope

Dino-Lite Digital Microscope The Industry Standard

	DINO-LITE 311S/313/411/413/423 SERIES	DINO-LITE 2XXX/3XXX SERIES	DINO-LITE 3613/3713/3715 SERIES	DINO-LITE 4013/4113/4023 SERIES	DINO-LITE AM 4115/4514*/4815* *SERIES	DINO-LITE WF/AF 4115*/4515*/4915** SERIES	DINO-LITE 7013/7023 SERIES	DINO-LITE 7115*/7515*/7815**/ 7915**/7025X** SERIES	DINO-LITE 73115*/73515*/ 73915** SERIES
Mac OS X 10.5 (Leopard)	Version 1.15 + Driver	Version 1.15	No	Version 1.15	No	No	Version 1.15	No	No
Mac OS X 10.6 (Snow Leopard)	Version 1.15 + Driver	Version 1.15	No	Version 1.15	No	No	Version 1.15	No	No
Mac OS X 10.7 (Lion)	Version 1.15 + Driver	Version 1.15	No	Version 1.15	No	No	Version 1.15	No	No
Mac OS X 10.8 (Mountain Lion)	Version 1.15 + Driver	Version 1.15	No	Version 1.15	No	No	Version 1.15	No	No
Mac OS X 10.9 (Mavericks)	Version 1.15 + Driver	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.10 (Yosemite)	Version 1.15 + Driver	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.11 (El Capitan)	No	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	No	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.12 (Sierra)	No	Version 1.24	No	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24
Mac OS X 10.13 (High Sierra)	No	Version 1.24	No	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24
Mac OS X 10.14 (Mojave)	No	Version 1.24	No	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24

Dino-Lite Digital Microscope Medical

	DINO-LITE MEDICAL DERMASCOPE MEDL4DW, MEDL7DW	DINO-LITE MEDICAL DERMASCOPE MEDL4DM, MEDL7DM	DINO-LITE MEDICAL EARSCOPE MEDL3E, MEDL4E, MEDL4EP	DINO-LITE MEDICAL TRICHOSCOPE MEDL3H, MEDL4HM/A/A5, MEDL4HVW, MEDL7HM/A/4	DINO-LITE MEDICAL CAPILLARYSCOPE MEDL4N PRO*, MEDL4N5 PRO*	DINO-LITE MEDICAL CAPILLARYSCOPE MEDL4N5	DINO-LITE MEDICAL IRISCOPE MEDL4R	DINO-LITE MEDICAL PODOSCOPE MEDL4FW
Mac OS X 10.5 (Leopard)	No	Version 1.15	Version 1.15	Version 1.15	No	Version 1.15	Version 1.15	No
Mac OS X 10.6 (Snow Leopard)	No	Version 1.15	Version 1.15	Version 1.15	No	Version 1.15	Version 1.15	No
Mac OS X 10.7 (Lion)	No	Version 1.15	Version 1.15	Version 1.15	No	Version 1.15	Version 1.15	No
Mac OS X 10.8 (Mountain Lion)	No	Version 1.15	Version 1.15	Version 1.15	No	Version 1.15	Version 1.15	No
Mac OS X 10.9 (Mavericks)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.10 (Yosemite)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.11 (El Capitan)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.12 (Sierra)	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24
Mac OS X 10.13 (High Sierra)	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24
Mac OS X 10.14 (Mojave)	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24	Version 1.24

* FLC and AMR compatible from version 1.22 and up

** FLC and AMR compatible from version 1.22. No EDOF/EDR functions available